

# Esthetic Tissue Formers

Esthetic Tissue Formers are made in a range of anatomical shapes for all positions and enables simple creation of the optimal emergence profile. The natural profile developed during provisionalization is matched perfectly in permanent restorative components. The design allows easy, fast and esthetic provisionals.

## Features & benefits

**Intended use** – Cement or screw retained provisional restorations and are recommended for single unit. Optionally, Esthetic Tissue Formers can be used as healing abutments.

**Packaging** – Sterile supplied with a titanium screw.

**Material** – Titanium and bondable polymer.

**Esthetiline Solution** – Esthetic Tissue Formers are part of the Esthetiline solution. The optimal result is achieved by choosing the same type of provisional and permanent restoration and the same position as during healing.

**Abutment modification** – Esthetic Tissue Formers may be easily modified by grinding with a bur chairside, intra- or extra-orally, or by the technician on a laboratory model.

**Indexing** – It is recommended that the implant is indexed prior to placing the Esthetic Tissue Former to ensure optimal placement and minimal preparation. All Esthetic Tissue Formers have a plane in the buccal direction and are non-rotational.

**Procedure** – Placement of Esthetic Tissue Formers at implant placement or abutment connection, guides the soft tissue and enables simple creation of the optimal emergence profile. For more information please refer to the Neoss Implant System Guidelines.

*Note: The trans-gingival section on Esthetic Tissue Formers is slightly smaller buccally than matching restorative components in order to provide additional soft tissue volume.*



## Assortment

Art. No.	Description
31340	Esthetic Tissue Former - Wide Incisor
31341	Esthetic Tissue Former - Narrow Incisor
31342	Esthetic Tissue Former - Canine
31343	Esthetic Tissue Former - Pre-molar
31344	Esthetic Tissue Former - Molar

