Glove Specifications & Packaging with Health Canada Covid-19 Interim Order Approval



Coronavirus disease (COVID-19)

Authorized medical devices for uses related to COVID-19: List of authorized medical devices other than testing devices

Device Name   ↑   •	Manufacturer 🚹 🖶	Class ↑↓	Technology 1  ↑  ↓	Date authorized ↑↓
Basic Vinyl Synmax Exam Gloves	Shandong Intco Medical Products Co. Ltd.	II	PPE	2020-05-28



Medical Devices Directorate Direction des instruments médicaux

Components/Parts/Accessories/Devices for this Licence
Les composantes, parties, accessoires et instruments médicaux pour cette homologation

#### BASIC VINYL SYNMAX EXAM GLOVES

Device ID/No de l'instrument: 1022102 Device Identifier / Identificateur de l'instrument (Model/Catalog Detail/No de modèle/Catalogue): BMPF3001 BMPF3002 BMPF3003



### Shandong Intco Medical Products Co., Ltd.

Add:No.9888 Qiwang Road,Naoshan Industrial Park
Qingzhou,Shandong,China
Web:www.intcomedical.com

### Vinyl Disposable Gloves Product Specification

Туре	Dispos	isposable Blue Synmax Vinyl Gloves,Powder Free											
Grade		Medi	cal	Surface	urface Smooth				Glove Length(mm)≥230				
Material	Naterial Poly-Vinyl Cloride		+ NBR LATEX	Packing		100pcs/box,10boxes/ctn Weight(g) 5.0±0.3(Medium Size)				100pcs/box,10boxes/ctn			
Cuff	Beaded		Shelf-life	е	5 years			Powder Level	≤2mg/glove				
				Physica	l Dimens	ion			Physical Property				
Product	Color	Size	Weight (g)	ILength	Palm Width (mm)	Cuff Thickness (mm)		Finger Thickness (mm)	Tensile Strength(MPa)	Elongation(%) Water Test AQL		Appearance AQL	
Disposable		S	4.5±0.3g	≥230	85±5	0.07±0.03	0.09±0.03	0.11±0.03		≥350	AQL2.5	Critical AQL1.0 Major AQL2.5 Minor 4.0	
Blue Synmax Vinvl	Blue	М	5.0±0.3g	≥230	95±5	0.07±0.03	0.09±0.03	0.11±0.03	>12		AQL2.5	Critical AQL1.0 Major AQL2.5 Minor 4.0	
Gloves, Powder		L	5.5±0.3g	≥230	105±5	0.07±0.03	0.09±0.03	0.11±0.03	≥13	≈35U	AQL2.5	Critical AQL1.0 Major AQL2.5 Minor 4.0	
Free		XL	6.0±0.3q	≥230	115±5	0.07±0.03	0.09±0.03	0.11±0.03			AQL2.5	Critical AQL1.0 Major AQL2.5 Minor 4.0	

We has established the Quality management system in according with ISO9001 and ISO 13485 standards

The certification which we have is 510K ,CE certification ,ISO9001 ISO13485,

Sampling Procedures according to the standard of ISO,ASTM D5250 ,Inspection Level for pysical dimensions and force /elongation at break is according to S-2 ,water test according to G-1





## Health Canada Covid-19 Interim Order Approval



Medical Devices Directorate Direction des instruments médicaux

### COVID-19 Medical Device Authorization for Importation or Sale

Authorization Reference Number: 315664 Numéro de référence de l'autorisation

 Issue Date:
 2020-05-28
 Date de délivrance:

 Amendment Date:
 2021-01-21,
 Date de modification:

 Reason for Amendment
 Correction to Device Class
 Raison de la modification

#### Device Class/Classe de l'instrument : 2

Pursuant to section 5 of the Interim Order Respecting the Importation and Sale of Medical Devices for Use in Relation to COVID-19, made by the Minister of Health on March 18, 2020, the medical device listed below is now authorized for sale or importation in Canada.

Each shipment of a COVID-19 medical device that is imported into Canada must be accompanied by a copy of this authorization document. Please ensure to highlight the Authorization reference number during the import declaration process to facilitate port entry without any delays.

This authorization is only valid for so long as the Interim Order Respecting the Importation and Sale of Medical Devices for Use in Relation to COVID-19 is in effect. Conformément à l'article 5 de l'Arrêté d'urgence concernant l'importation et la vente d'instruments médicaux relatifs au Covid-19, réalisé par la ministre de la Santé le 18 mars 2020, les instruments indiqués cidessous sont présentement autorisés pour la mise en vente ou l'importation au Canada.

Autorisation d'importation ou de

mise en vente d'un instrument

médical relatif au COVID-19

Tout envoi d'un instrument médical relatif au COVID-19 doit être accompagné d'une copie de la présente autorisation. Veuillez vous assurez de souligner le numéro de référence de l'autorisation durant le processus de déclaration d'importation pour faciliter l'entrée sans delais aux points de contrôle frontelier.

Cette autorisation est uniquement valide tant que l'Arrété d'urgence concernant l'importation et la vente d'instruments médicaux relatifs au Covid-19 est en vieueur. ou l'autorisation est annulée.

139268


#### Device Name(s) Nom de l'instrument

#### BASIC VINYL SYNMAX EXAM GLOVES

Name & Address of Authorization Holder/Nom & adresse du titulaire de l'autorisation

SHANDONG INTCO MEDICAL PRODUCTS CO., LTD. QIWANG ROAD, NAOSHAN INDUSTRIAL PARK, QINGZHOU WEIFANG, SHANDONG CHINA

Application Number: Manufacturer ID: Numéro de la demande: 315664 Identificateur du fabricant:



Application Number:

Numéro de la demande



Medical Devices Directorate Direction des instruments médicaux

Components/Parts/Accessories/Devices for this Licence
Les composantes, parties, accessoires et instruments médicaux pour cette homologation

#### BASIC VINYL SYNMAX EXAM GLOVES

Manufacturer ID:

Identificateur du fabricant:

139268

Device ID/No de l'instrument: 1022102 Device Identifier / Identificateur de l'instrument (Model/Catalog Detail/No de modéle/Catalogue): BMPF3001 BMPF3002 BMPF3003 BMPF3004

Manufacturer ID: 139268 Application Number: 315664

Identificateur du fabricant: Numéro de la demande:

### ASTM D5250-19 Test Report



Test Report

No.: QDHL2004002960MD-01

THIS REPORT IS TO SUPERSEDE TEST REPORT NO.: QDHL2004002960MD, DATE: APR.29,2020 THE ORIGINAL REPORT SHALL BECOME INVALID AS OF THE DATE OF ISSUANCE OF THIS

Client name SHANDONG INTCO MEDICAL PRODUCTS CO.,LTD Client address

NO.9888, QIWANG ROAD, NAOSHAN INDUSTRY PARK, QINGZHOU, SHANDONG, CHINA

SYNMAX VINYL EXAM GLOVES\* Sample Description

Lot No. : NOT PROVIDED Lot Size : NOT PROVIDED

As above test item and its relevant information regarding to the submission are provided and confirmed by the applicant. SGS is not liable to either the test item or its relevant information, in terms of the accuracy,

suitability, reliability or/and integrity accordingly.

Sample Receiving Date : APR.14,2020 Test Performing Date APR.14,2020 TO APR.29,2020

Test Performed SELECTED TEST(S) AS REQUESTED BY APPLICANT

ASTM D 5250-19 STANDARD SPECIFICATION FOR POLY

(VINYL CHLORIDE) GLOVES FOR MEDICAL APPLICATION

(CLAUSE 6.1.2, 6.1.3, 6.1.4, 6.1.5)

Test Result(s) PLEASE REFER TO THE FOLLOWING PAGE(S)

Conclusion : THE SUBMITTED SAMPLE MET THE TEST REQUIREMENT.

Remark: Unless otherwise stated the results shown in this test report refer only to the sample(s) tested.

SGS-CSTC Standards Technical Services (Qingdao)

Co., Ltd.

Test Requested

Jessica Gao Approved Signatory



SCR Contor No. 123 Thurston Road Landbur District Closelan China 26010

Member of the SGS Group (SGS SA)



No.: QDHL2004002960MD-01

Date: SEP.22.2020 Page: 2 of 5

Test Conducted

ASTM D 5250-19 Standard Specification for Poly (vinyl chloride) Gloves for Medical Application

Number of t	est sample		: 231 Piece(s)		
Clause	Test Items			Result	Note
6.1.2	Freedom fr	rom holes		Pass	#1
6.1.3	Physical di	mensions		Pass	#2
6.1.4	Physical pr	operty charac	cteristics	Pass	#3
6.1.5	Powder res	sidue for pow	der free gloves	Pass	#4
Notes	: #1 S	iee result 1			
	#2 0	on specific 2			

See result 2

See result 3 The average mass of powder per glove is 0.04mg

Test Result:

Result 1: Freedom from Holes

Sample Quantity: 200 Ac: 10 Re: 11 Found: 0



1909 Contar No. 163 Zharthou Road Landbou District Claration China 20010

Member of the SGS Group (SGS SA)

## ASTM D5250-19 Test Report



No.: QDHL2004002960MD-01 Date: SEP.22,2020 Page: 3 of 5

No.: QDHL2004002960MD-01

No.: QDHL2004002960MD-01

Date: SEP.22,2020 Page: 5 of 5

### Force at break after aging

	Size: M					
Sample No.	Tensile strength (Mpa)	Ultimate Elongation (%)				
1	19	370				
2	20	409				
3	18	382				
4	19	374				
5	20	378				
6	18	356				
7	16	347				
8	19	360				
9	18	372				
10	16	307				
11	18	368				
12	17	340				
13	18	357				
Standard requirement	≥11	≥300				
Found	0	0				

- 1. The declaration of conformity is only based on the actual value of laboratory activity, measurement uncertainty of the results not take into account.
- The marked ▲ part is the modification information.

#### Sample Photo:



SGS authenticate the photo on original report only

""End of Report""



SSS Contor, No. 143, Zhuzhou Road, Lanshen District, Gingdon, China 26610

Member of the SGS Group (SGS SA)

Result 2: Physical dimensions

Sample Quantity: 13 Ac: 1 Re: 2 Found: 0

	Size: M									
			Median value/mm							
Sample No.	Length/mm	Width/mm	Thickness-finger	Thickness-palm						
1	238	96	0.082	0.091						
2	237	98	0.094	0.098						
3	237	98	0.092	0.092						
4	240	96	0.093	0.097						
5	232	97	0.092	0.096						
6	236	98	0.092	0.097						
7	230	96	0.094	0.093						
8	234	96	0.083	0.095						
9	235	97	0.081	0.094						
10	239	98	0.113	0.093						
11	240	97	0.092	0.095						
12	242	98	0.135	0.091						
13	239	98	0.143	0.095						
Standard requirement	≥230	95±5	≥0.08	≥0.08						
Found	0	0	0	0						



Result 3: Physical property characteristics

Sample Quantity: 13 Ac: 1 Re: 2 Found: 0

Force at break before aging

	Size: M				
Sample No.	Tensile strength (Mpa)	Ultimate Elongation (%)			
1	16	313			
2	16	394			
3	17	344			
4	19	364			
5	18	365			
6	20	395			
7	19	372			
8	17	346			
9	21	404			
10	18	349			
11	18	365			
12	19	366			
13	19	368			
Standard equirement	≥11	≥300			
Found	0	0			



SGS Center, No. 143, Zhuzhou Road, Laeshen Diebirt, Gingdon, China 26610

SGS Contac No. 143, Zhuzhou Road, Lanshen District, Gington, China 26810

Member of the SGS Group (SGS SA)

Member of the SGS Group (SGS SA)

### ASTM D6499 Test Report for The Immunological Measurement of Antigenic Protein



Testing, Development, Problem Solving,

October 20, 2020

\*TEST REPORT\*

PN 156280 Check

CHEMICAL ANALYTICAL SERVICES

Prepared For:

Cecily Sheng Intco Medical Industries Inc. 805 Barrington Avenue Ontario CA 91764

Prepared By

Erin Samples, Sr. Technician, Microbiology

Chemical/Analytical Services

Approved By:

Ana Barbur, Vice President Chemical/Analytical Services

Rev 110119

An A2LA ISO 17025 Accredited Testing Laboratory - Certificate Numbers 255.01 & 256.02 ISO 9001-2015 Registered

ISO 9001:2015

This report is for the exclusive use of the client(s) to which it is addressed and shall not be reproduced, except in full, without written permission of Akron Rubber Development Laboratory. Inc. (ARDL). The information provided herein applies to the specific material, products, or processes tented or evaluated. Statements of Conformity instituted herein are based on the option as agreed in the Client testing contract. No warranty of any kind is been construed or impired. The liability of ARDL is the limited to the mission of coefficient or processes and the conformation of the liability of ARDL is the limited to the limited on the information conformation.

> www.ardl.com | 2887 Gilchrist Rd. | Akron, Ohio 44305 | answers@ardl.com | Toll Free (800) 830-ARDL. Fax (330) 794-6610 | Worldwide (330) 794-6600



Testing, Development, Problem Solving,

October 20, 2020 Cecily Sheng Intco Medical Industries Inc.

Page 2 of 4 PN 156280

Analytical testing on samples submitted by the above referenced customer

One (1) glove sample identified as Blue Powder Free Synmax Size M Gloves

**Decision Rule 1** 

ASSAY PROCEDURE

Sample Extractions. The sample was weighed, measured, and cut prior to extraction within a polypropylene extraction vessel. The extraction buffer used was 50 mM phosphate pH 7.4 at a ratio 5 ml of buffer per gram of sample. The extraction was carried out at room temperature for two hours with agitation. The sample was removed and the extract centrifuged at 500 xg for 15 minutes to pellet particulates. The cleared extract was then used in the assay.

ASTM D6499 ELISA Inhibition Assay. The standard and test samples are serially diluted in a 96 well plate, after which an equal volume of diluted rabbit anti-latex polyclonal antibody is added and the plate incubated for 2h at 37°C. One hundred microliters of sample from each well is transferred to the corresponding well of a plate coated with Hevea NRL and blocked with non-fat dry milk and incubated for 2h at 37°C. The plates are then washed and a 100 µl solution of Goat anti-Rabbit IgG conjugated with the enzyme HorseRadish Peroxidase (HRP) is added and incubated for 1h at 37°C. Plates are washed and a 100-ul solution of the substrate OPD is added to each well and color allowed to develop. The reaction is stopped by the addition of 50 ul of 4N H<sub>2</sub>SO<sub>4</sub>. The plate is then read at 490 nm. Protein values are determined by interpolation from a

RESULTS: The sample identified as Blue Powder Free Synmax Size M Gloves tested below the detection limit of the assay at ≤0.2 µg/gm and ≤0.1 µg/dm2.

#### **ASTM D6499 Test Certificate**

Sample Description	Weight (gm)	Area (dm²)	Extract Vol. (ml)	Assay Conc. (µg/ml)	Antigenic Protein (µg/gm)		Antigenic Protein (ug/dm²)	
Blue Powder Free Synmax Size M Gloves	4.9	8.7	24.5	b.d.	<	0.2	<	0.1

Where b.d.= below detection, (0.03 µg/ml)

\*Note: The results given in this report relate only to the items tested. This report cannot be reproduced except in full without the written consent of ARDL

\*ARDL is ISO 17025 accredited by A2LA for the test methods listed on the certificates referenced on page one. Unless specified, the current specification version is used. NOTE: Non-ISO 17025 accredited test methods are designated with the ^ symbol to differentiate from ISO 17025 accredited methods in the body of the test report.

> www.ardl.com | 2887 Gilchrist Rd. | Akron. Ohio 44305 | answers@ardl.com | Toll Free (800) 830-ARDL Fax (330) 794-6610 | Worldwide (330) 794-6600

## ASTM D6499 Test Report for the Immunological Measurement of Antigenic Protein

October 20, 2020 Cecily Sheng Intco Medical Industries Inc.



Sample - Blue Glove Powder Free Synmax Size M

Page 3 of 4 PN 156280

October 20, 2020 Cecily Sheng Intco Medical Industries Inc.

Page 4 of 4 PN 156280

#### **Decision Rules**

Rule 1. This is the way test results have traditionally been reported by ARDL. If ARDL runs a test for you that has pass/fail requirements, ARDL will report the values observed and then state "Pass" or "Fail", based on those values only. By default, ARDL will apply this rule to all Category I tests and those tests which are not on ARDL's Scope of Accreditation.

Rule 2. This rule takes into account the calculated measurement uncertainty of test results generated. Every test and piece of test equipment has an inherent amount of measurement uncertainty associated with it. Rule 2 establishes "Guard Bands" where the measurement uncertainty value is added to the Minimum Passing requirement and is subtracted from the Maximum Passing requirement. The Pass/Fail requirements thus become tighter and customers may be more "Certain" of

Rule 3. This rule also takes into account measurement uncertainty but does not set up guard bands. Rule 3 may be used when values are reported, but there is no Pass/Fail requirement called out in the test specification. Rule 3 simply states that the measurement uncertainty is reported to the customer, along with the testing result generated, and the customer decides if the results are suitable for their purposes.

#### Report Revision Log

Report Revision 10-20-20

Description

Chemical/Analytical Services

Erin Samples, Sr. Technician, Microbiology

Ana Barbur, Vice Presiden Chemical/Analytical Services

ES/AB/tkr

## ASTM D1671 Test Report for Penetration of Phi-X174 Bacteriophage



Note: (1) General Terms & Conditions as mentioned overleaf. (2) The results relate only to the items tested.(3) The test report shall not be reproduced except in full without the written approval of the jaboratory.(4) Without the agreement of the jaboratory , the client is no authorized to use the test results for unapproved propaganda.

TÜV SÜD Produots Testing (Shanghai) Fex: +86 (21) 6037 6345

Phone: +86 (21) 6037 6375 Fmail: food chem@tuv-sud.cn

Regional Head Office: Jiangsu TÜV Products Service Co. Ltd. Shannhai Bransh 200 070 P.R.China

Test Report No.: 721659719-1 Report Date: 26 November 2020



#### RECEIPT DATE / TEST DATE

16-Nov-2020/ 16-Nov-2020

#### THE FOLLOWING SAMPLE(S) WAS/WERE SUBMITTED

#### BY/ ON BEHALF OF THE CLIENTS AS:

Sample Name: Synmax Vinyl Gloves

Sample Type: Sample Batch:

Shandong Intco Medical Products Co., Ltd. Manufactory:

SAMPLE NO.	DESCRIPTION	PHOTOGRAPH
721659719	Blue disposable gloves	(ALLEGAL)

#### TEST METHOD(S)

Penetration of Phi-X174 Bacteriophage Test

- with reference to ASTM F1871/F1871M-2013 Standard Test Method for Resistance of Materials Used in Protective Clothing to Penetration by Blood-Borne Pathogens Using Phi-X174 Bacteriophage Penetration as a Test System

#### REQUIREMENT

- Exposure Procedure: B

Retaining Screen:Stainless steel mesh screen Diameter of the steel wire: 0.80mm Retaining Screen pore: 2.0mm\*2.0mm Sampling Size: 75mm×75mm

Negative control: Polyethylene material

Positive control: 0.04 µm microporous membrane

Prior to testing, condition all test specimens and controls for a minimum of 24 hours at (21±5)°C and 30%~80% relative humidity.

#### TEST ORGANISM(S)

Bacteriophage ATCC 13706-B1

#### PROCEDURE

- 1.1. Test three specimens representing each material type to be tested.
- 1.2. With the sterile cell placed horizontally on the laboratry bench, insert the specimen in the test cell with the back of the hand and the palm of the hand toward the cell reservoir.
- 1.3. Assemble the test cell. Torque the bolts in the test cell to 13.6 N·m each.
- 1.4. With the cell remaining in the horizontal position on the laboratory bench, place a 2.0 µL aliquot of the Phi-X174 bacteriophage in bacteriophage nutrient broth, containing a total of 900-1200 PFU, near the middle of each piece of test specimen.

B-3/4, No.1999 Du Hul Road, Minhang

Email: food chemilitus-sud co

Regional Head Office: Jiangsu TÜV Products Bervice Co. Ltd. Shannhai Branch 200 070 P.R.China

Page 2 of 4



## ASTM D1671 Test Report for Penetration of Phi-X174 Bacteriophage

Test Report No.: 721659719-1 Report Date: 26 November 2020



- 1.5. Prepare a control by adding a 2.0 µL aliquot from the same suspension directly into 5.0 mL of sterile bacteriophage nutrient broth.
- 1.6. After 60 min, quantitatively assay by adding 5.0mL of sterile bacteriophage nutrient broth onto the surface of the specimen.
- 1.7. Calculate the ratio of the control assay titer to the test material assay titer using the following equation:

control assay titer (PFU/mL)

ratio=test material assay titer (PFU/mL) =1.2

- 1.8. Record the initial titer of the Phi-X174 bacteriophage challenge suspension used for the test. (2±1)x10<sup>8</sup> PFU/mL times the ratio calculated.)
- Test procedure
- 2.1. Prepare the bacteriophage challenge suspension (40-44 mN/m) for the test.
- 2.2. With the sterile cell placed horizontally on the laboratry bench, insert the specimen in the test cell with the back of the hand and the palm of the hand toward the cell reservoir.
- 2.3. Assemble the test cell. Torque the bolts in the test cell to 13.6 N·m each.
- 2.4. Mount the test cell in the test apparatus in a vertical position and close the drain valve.
- Penetration test: If liquid appears to penetrate through the test specimens at anytime during the test, terminate the test.
  - (1) Carefully fill the test cell reservoir with approximately 60 mL of the Phi-XI74 bacteriophage challenge suspension
  - (2) Step1: Observe for 5 min at 0 psi.
  - Step2: Slowly increase the pressure to 2.0 psi at the rate of no more than 0.5 psi/s, keep the pressure at 2.0 psi, observe for 1 min.
  - Step3: Slowly decrease the pressure to 0 psi at the rate of no more than 0.5 psi/s, observe for 54 min.
- (3) At the end of the time period, open the drain valve and drain the test cell of the bacteriophage challenge suspension. Dilute and assay the bacteriophage concentration. 2.6. Specimen surface assay procedure
  - (1) With the sterile cell placed horizontally on the laboratory bench. Slowly add 5.0 mL of sterile bacteriophage nutrient broth onto the exposed surface of the specimen.
  - (2) Gently swirl the test cell for approximately 1 min. Assay the liquid soon after collecting.
- 2.7.Remove the specimen from the test cell and prepare the test cell for sterilization.
- Test controls
  - 3.1. The negative control was negative for bacteriophage penetration.
  - The positive control was positive for bacteriophage penetration.
  - 3.3. Record the final titer of the bacteriophage challenge at the conclusion of the 60 min test.

Test Report No.: 721659719-1 Report Date: 26 November 2020



#### TEST RESULT(S)

		Initial titer	Final titer	Test Results					
Test Items		PFU/mL	PFU/mL	Step1	Step2	Step3	Assay titer (PFU/mL)	PASS/FAIL	
Penetration of Phi-X174 Bacteriophage	Control(+)	1.8×10 <sup>8</sup>	1.7×10 <sup>8</sup>	None Seen	Seen			Acceptable	
	Control(-)	1.8×10 <sup>8</sup>	1.7×10 <sup>8</sup>	None Seen	None Seen	None Seen	<1	Acceptable	
	-1	1.8×10 <sup>8</sup>	1.7×10 <sup>8</sup>	None Seen	None Seen	None Seen	<1	PASS	
	-2	1.8×10 <sup>8</sup>	1.7×10 <sup>8</sup>	None Seen	None Seen	None Seen	<1	PASS	
	-3	1.8×10 <sup>8</sup>	1.7×10 <sup>8</sup>	None Seen	None Seen	None Seen	<1	PASS	

Note:

PFU: Plaque Forming Unit.

2. This report is for internal use only such as internal scientific research ,education, quality control, product







ChemicalMicrobiology Laboratory: TUV 8UD Products Testing (Shanghai) Co., Ltd B-344, No.1999 Du Hul Roed, Minhang District Shanghel Phone : +86 (21) 6037 6375
Fax : +86 (21) 6037 6345
Email: food.chem@tuv-sud.cn
Webpage: www.huv-sud.cn

Regional Head Office: Jiangsu TÜV Produots Service Co.,Ltd. Shanghai Branch. No.151 Heng Tong Rood Shenghel 20070 P.R.Chine

ChemicalMicrobiology Laboratory: TÜV SÜD Products Testing (Shanghai Co., Ltd B-3/4, No.1999 Du Hul Roed, Minheng District Phone: +86 (21) 6037 6375

Fex: +86 (21) 6037 6345

Regional Head Office: Jiangsu TUV Products Service Co.,Ltd. Shanghai Branch. No.151 Heng Tong Roed Shanghal 200 070 \_P.R.Chine

Page 4 of 4

Glove Specifications & Packaging with Health Canada Covid-19 Interim Order Approval

## Applications



Medical Purpose / Examination



Industrial purpose / PPE



Laboratory



Healthcare and nursing



General housekeeping



IT Industry






















<sup>\*</sup> The last "\*" means Size:4-S 5-M 6-L 7-XL.

ISO 9001:2015, ISO 13485:2016 and MDSAP Audit Certification



## CERTIFICATE OF REGISTRATION

This is to certify that the management system of:

### **Shandong Intco Medical** Products Co., Ltd.

Qiwang Road, Naoshan Industrial Park, Qingzhou City, Shandong Province,

has been registered by Intertek as conforming to the requirements of:

### ISO 9001:2015

The management system is applicable to:

Manufacture of disposable non sterile NBR (Nitrile Butadiene Rubber) and PVC (Poly Vinyl Chloride) gloves

Unified Social Credit Identifier: 91370781561439654L

#### Certificate Number: 111812005

Initial Certification Date: 27 January 2013

Date of Certification Decision: 03 January 2019

### Issuing Date:

03 January 2019

#### Valid Until: 27 January 2022



#### Calin Moldovean

President, Business Assurance

Intertek Certification Limited, 10A Victory Park, Victory Road, Derby DE24 8ZF, United

Intertek Certification Limited is a UKAS accredited body under schedule of accreditation no. 014.



## CERTIFICATE OF REGISTRATION

This is to certify that the management system of:

### **Shandong Intco Medical** Products Co., Ltd.

Main Site: Qiwang Road, Naoshan Industrial Park, Qingzhou City, Shandong Province, 262500 P.R. China

has been registered by Intertek as conforming to the requirements of:

### ISO 13485:2016

The management system is applicable to:

Manufacture of non-sterile NBR (Nitrile Butadiene Rubber) and PVC (Poly Vinyl Chloride) patient examination gloves.

Certificate Number: 0086238

Initial Certification Date 28 April 2014

Certificate Issue Date: 17 January 2019

Certificate Expiry Date: 31 December 2021





Calin Moldovean President, Business Assurance

Intertek Testing Services NA Ltd., 1829, 32nd avenue, Lachine, QC, HST 3J1,



# **CERTIFICATE OF REGISTRATION**

This is to certify that the management system of:

## **Shandong Intco Medical** Products Co., Ltd.

(DUNS # 42-128-5853)

Main Site: Qiwang Road, Naoshan Industrial Park, Qingzhou City, Shandong Province, 262500 P.R. China

has been registered by Intertek, an MDSAP recognized auditing organization, as conforming to the requirements of:

### ISO 13485:2016

Canada: Medical Devices Regulations - Part 1- SOR 98/282

United States: 21 CFR 820, 21 CFR 803, 21 CFR 806, 21 CFR 807 (Subparts A to D)

#### The management system is applicable to:

Manufacture of non-sterile NBR (Nitrile Butadiene Rubber) and PVC (Poly Vinyl Chloride) patient examination gloves.

Certificate Number

Initial Certification Date:

2019-01-25


Certification Effective Date

2019-01-25

0086239

Certification Expiry Date: 2022-01-24






President, Business Assurance

Intertek Testing Services NA, Inc. 900 Chelmsford Street Lowell, MA, USA 01851














the inverse of this carrificate, intertels assumes no Eablity to any party other than to the Client, and then only in accordance with the agreed upon Certificati

Genuine access to Intco via official communication channels



Our official email address suffix ends with @intco.com. It's the only email address we have.

(office) (cell)

@intco.com

# Intco Vinyl Synthetic Exam Gloves with Health Canada Approval

Sample sent directly from Intco Factory



