



# 3D High Content Imaging of Optically-Cleared Microtissues for Screening Islet Cell Proliferation

Thomas S. Villani, Graeme Gardner, Michael Johnson, Burcak Yesildag, Joan Mir  
Visikol Inc, 295 US Highway 22 E, Whitehouse Station, NJ, 08889 USA  
InSphero AG, Wagistrasse 27, Schlieren, 8952 Switzerland

## Highlights

When combined with high content confocal microscopy, Visikol® HISTO-M™ allows the interior environment of InSphero 3D InSight™ Microtissues to be imaged, dramatically improving quantification of proliferating cells in pancreatic islets.

- Standardized 3D InSight Human Islet Microtissues recapitulate native architecture, function and proliferative capacity
- High-Content Imaging along with tissue clearing enables full quantification of specific cell populations in 3D cell culture

## Introduction

The use of *in vitro* three-dimensional (3D) cell cultures has increased dramatically for drug discovery because 3D cell culture models more accurately mimic the *in vivo* environment compared to traditional monolayer cultures<sup>1</sup>. However, current imaging-based analysis of these 3D cultures predominantly relies upon techniques developed originally for 2D cell culture. Due to the thickness of 3D cell cultures, typically >100 µm, light scattering does not permit imaging of the center of the spheroid<sup>2</sup>. This technical limitation introduces a sampling bias in imaging analysis, as only the exterior cells can be imaged. This problem must be solved to obtain the most useful and accurate survey of the cellular environment and response of the spheroid.

Understanding pancreatic β-cell function in the regulation of glucose homeostasis and metabolism is of paramount importance in the research of the most common human metabolic disease – diabetes<sup>3</sup>. Although isolated primary islets are considered the gold standard tool of diabetes research, their experimental use is

limited due to their inherent heterogeneity in size, cellular composition and purity, as well as their rapid decline in functionality and viability *ex vivo*.

The ability of the small molecule harmine to activate β-cell proliferation *in vivo* is a well-established phenomenon and, for this reason, has initiated interest developing analog compounds to treat diabetes<sup>4</sup>. This in turn has given rise to interest in the development of high-throughput screening (HTS) assays to survey the efficacy of these drugs in *in vitro* models of pancreatic islets.

In this study, we combine the 3D InSight™ Human Islet Microtissues provided by InSphero and the Visikol HISTO-M tissue clearing reagent to illustrate the power of tissue clearing on high content screening (HCS) of novel drugs targeting β-cell proliferation.

## Materials and Methods

### Reagents and cell culture

3D InSight™ Human islet microtissues were obtained from InSphero Inc. Antibodies and fluorescent dyes and

## Application Note

antibodies were obtained from Invitrogen (Thermo Scientific). Harmine was obtained from Sigma Aldrich.

3D InSight™ Islet Microtissues were prepared from human islets (Age 29, Male, Hispanic, BMI 22.17, HbA1c 5.5%) and delivered in an Akura™ 96 well plate. Medium exchanges were performed every 2 days to maintain the microtissues in culture. On day 19, a 4.5-day treatment with 10  $\mu$ M Harmine was started, with re-dosing after 2.5 days. Harmine was dissolved in DMSO to make a 10 mM stock solution, and a vehicle control was treated alongside for comparison.

### Fixation and immunolabeling

On day 23, microtissues were fixed using 4% PFA, followed by washing in PBS to remove fixative. Microtissues were treated with methanol, followed by 20% DMSO/methanol, followed by 0.2% Triton X-100 in PBS to improve penetration of antibodies. Microtissues were blocked with 10% donkey serum, then incubated with guinea pig anti-PDX1 (1:150 dilution, Abcam Cat# ab47308), rabbit anti-Ki67 antibody (1:200 dilution, Abcam Cat# ab15580), and with goat anti-rabbit AlexaFluor568 and goat anti-guinea pig AlexaFluor488 conjugated secondaries (1:250 dilution, available from Invitrogen) to label  $\beta$ -cells and proliferating cells. Nuclei were counterstained with DAPI (Molecular Probes).

### Clearing and high throughput imaging of microtissues


Microtissues were dehydrated with methanol, then cleared with Visikol HISTO-M. Imaging of spheroids was conducted in InSphero Akura™ 96 well plates using a CX7-LZR High Content Screening platform confocal imager (ThermoFisher). Multi-channel z-stacks were collected for each tissue, using 5  $\mu$ m steps at 20X magnification. Images were processed using automated ImageJ macros, and cells were counted using CellProfiler.

## Results and Discussion

### Drug-induced $\beta$ -Cell Proliferation

The islet microtissues were treated with harmine at a concentration known to induce significant proliferation in dispersed islets cells (10  $\mu$ M). After intervention with harmine, the effects on islet function were studied. Harmine did not alter the basal levels of insulin secretion at 2.8 mM glucose. Similarly, stimulated insulin secretion

levels were maintained in harmine treated islets, reflecting unaltered  $\beta$ -cell function. However, harmine induced a significant reduction in the levels of total insulin content (Figure 1).




**Figure 1.** Functional performance in islet spheroids. (Top) Total Insulin levels and (Bottom) secreted insulin in basal (2.8 mM) and stimulated (16.7 mM) glucose levels.

### Detection and Quantification of Proliferating $\beta$ -Cells

In order to probe the effect of target compounds on the total insulin-producing  $\beta$ -cell population within pancreatic islet microtissues, it is necessary to be able to quantify that sub-population of cells within the tissue and determine whether or not a significant increase in proliferation has occurred. For this, PDX1 (transcription factor associated with  $\beta$ -cell maturation) was identified as the appropriate biomarker label to detect  $\beta$ -cells and anti-Ki67 was used for labeling the proliferating sub-population along with DAPI to counterstain for all nuclei. The population of proliferating  $\beta$ -cells as well as total proliferating cells within each islet microtissue were then quantified by colocalization of the three signals. Both PDX1 and Ki67 are localized within the nucleus of the cell, making identification and colocalization with one


another along with the DAPI signal a straightforward analysis.

The image montages for both control and harmine-treated tissues (Figure 1) show the number of cells positive for Ki67 increases substantially between the two example image sets. A total of 12 tissues for each condition were surveyed such that statistical analysis could be performed and the summarized average total cell counts visualized (Figure 2). The tissues that were treated with harmine on average exhibit a 10-15X fold change increase in proliferating cells (with p-values for Ki67+ only =  $10^{-6}$ , and for Ki67+/PDX1+ =  $10^{-5}$ ).




**Figure 2.** Montages of z-stack images taken from CX7-LZR HCS platform showing individual cells can be counted from each section throughout the tissue, (DAPI = Blue, PDX1 = Green, Ki67 = Red). A) Vehicle control showing typical  $\beta$ -cell distribution in islet microtissue and minimal proliferating cells; B) Harmine (10  $\mu$ M) treated microtissue showing greatly increased proliferating cell population.

### Total Cell Counts



### Percent Proliferation



**Figure 3.** Summary of average cell totals counted for control and harmine-treated microtissues. (Top) Total cell counts are for all cells (blue), and for the  $\beta$ -cell sub-population (green). (Bottom) Percent proliferating cells for all (blue) and  $\beta$ -cell sub-population (green).

The average percentage of proliferating cells for the islet microtissues increases significantly as well, going from  $\sim 0.5\%$  proliferating  $\beta$ -cells in the control sample, to  $\sim 9.8\%$  proliferating  $\beta$ -cells in the harmine-treated samples, consistent with what has been observed previously<sup>4</sup>. The distribution of proliferating cells within each microtissue can be elicited from the high-content imaging data, showing that the proliferating population is dispersed relatively uniformly throughout the tissue, with cells on the interior having equal chance of being affected by harmine as those on the periphery.

When conducting confocal imaging of microtissues without clearing, only the outermost cells are detected due to light scattering that reduces signal to noise. Use of a tissue clearing agent greatly increases the number of cells detectable by high content imaging and eliminates

## Application Note

the need to do time-consuming and technically challenging IHC of microtissues. The non-adherent surface coating of InSphero assay plates was chemically compatible with Visikol HISTO-M. The clearing process takes only minutes and is done within the well plate.

## Conclusions

- Using Visikol HISTO-M coupled with confocal imaging, light scattering and out-of-focus light is reduced, allowing for comprehensive profiling of the interior of each islet microtissue.
- Visikol HISTO-M is compatible with immunolabeling and fluorescent staining of microtissues.
- Clearing microtissues with Visikol HISTO-M enables visualization of the interior of microtissues.
- Harmine induces a 10-15X fold change in the number of proliferating cells within 3D InSight™ Human Islet Microtissues, comparable to *in vivo* results.

## References

1. LaBarbera, D. V., Reid, B. G., & Yoo, B. H. (2012). The multicellular tumor spheroid model for high-throughput cancer drug discovery. *Expert opinion on drug discovery*, 7(9), 819-830.
2. Kabadi, P. K., Vantangoli, M. M., Rodd, A. L., Leary, E., Madnick, S. J., Morgan, J. R., & Boekelheide, K. (2015). Into the depths: Techniques for in vitro three-dimensional microtissue visualization. *Biotechniques*, 59(5), 279.
3. Bell, G.I.B., Polonsky, K.S. (2001). Diabetes mellitus and genetically programmed defects in beta-cell function. *Nature* 414, 788-791.
4. Wang, P., Alvarez-Perez, J-C., Felsenfeld, D.P., Liu, H., Sivendran, S., Bender, A., Kumar, A., Sanchez, R., Scott, D.K., Garcia-Ocana, A., Stewart, A.F. (2015). Induction of human pancreatic beta cell replication by inhibitors of dual specificity tyrosine regulated kinase. *Nat. Med.*, 21(4), 383-388.

## Getting Started

1. Email us at [info@visikol.com](mailto:info@visikol.com) to discuss your research project.
2. Visit [www.insphero.com](http://www.insphero.com) to learn more about pancreatic islet microtissues.
3. Collaborations begin with pilot projects which can then be scaled up according to your requirements.
4. We work end-to-end to help you take advantage of 3D High Content Imaging assays in your workflow.