



neoss®

NeoGen® Resorbable Membranes and dressings




Easy to use,
clinically proven

Intelligent Simplicity

Resorbable membranes

NeoGen® Collagen Flex

NeoGen® Collagen Flex Resorbable Membranes are a strong, conformable collagen barrier membrane manufactured from purified porcine peritoneum tissue. It is designed to provide a very drapable yet very strong membrane.




Product Features

- ✓ Resorbable in 3–4 months
 - ✓ High mechanical strength
 - ✓ Highly purified intact peritoneum tissue
-

Applications

- In augmentation around implants placed in immediate extraction sockets or delayed extraction sockets.
- Localized ridge augmentation for later implant placement.
- Alveolar ridge reconstruction for prosthetic treatment.
- Guided bone regeneration in dehiscence defects and guided tissue regeneration procedures in periodontal defects.



Scanning Electron Micrograph of the dry NeoGen® Collagen Flex (cross section) at 50x magnification

NeoGen® Collagen Firm

NeoGen® Collagen Firm Resorbable Membranes have optimized characteristics to provide periodontal and dental surgeons with the ideal balance of properties to effectively address a host of clinical indications and surgical procedures.

NeoGen® Collagen Firm is manufactured from porcine dermis tissue.




Product Features

- ✓ Resorbable in 6–9 months
- ✓ Highly purified intact dermis tissue
- ✓ Balanced characteristics of handling, strength and resorption time
- ✓ High mechanical strength

Applications

- In augmentation around implants placed in immediate extraction sockets or delayed extraction sockets.
- Localized ridge augmentation for later implant placement.
- Alveolar ridge reconstruction for prosthetic treatment.
- Guided bone regeneration in dehiscence defects and guided tissue regeneration procedures in periodontal defects.



Scanning Electron Micrograph of the dry NeoGen® Collagen Firm (cross section) at 50x magnification

NeoGen® Resorbable Collagen Wound Dressings

NeoGen® Resorbable Collagen Wound Dressings are absorbent, porous, collagen matrices engineered from purified collagen derived from bovine dermis tissue.

Essentially resorbs within 30 days of placement.



Product Features

- ✓ Applied directly to the wound and protects the wound and delicate new tissue
- ✓ Aids in wound healing
- ✓ Available in tape, foam and plug form

Indications

Oral ulcers (non-infected or viral)	✓
Burns	✓
Suture sites	✓
Extraction sites	✓
Surgical wounds	✓
Traumatic wounds	✓

Resorbable membranes and dressing assortment

	ARTICLE NO.	SIZE
NeoGen® Collagen Flex	64050	15 mm x 20 mm
	64051	20 mm x 30 mm
	64052	30 mm x 40 mm
NeoGen® Collagen Firm	64055	15 mm x 20 mm
	64056	20 mm x 30 mm
	64057	30 mm x 40 mm
NeoGen® Resorbable Collagen Wound dressings	64060	25 mm x 75 mm x 1 mm - 10pcs/pkg
	64061	20 mm x 40 mm x 3 mm - 10pcs/pkg
	64062	10 mm x 20 mm - 10pcs/pkg

CE 2797



Manufacturer: Collagen Matrix, Inc., 15 Thornton Road, Oakland, NJ 07436 USA
Phone: 201-405-1477 | Facsimile: 201-405-1355 | www.collagenmatrix.com

Authorized EU Representative: MDSS GmbH, Schiffgraben 41, 30175 Hannover, Germany

Distributed by: neoss.com

Neoss Ltd. Windsor House, Cornwall Road, Harrogate, HG1 2PW, United Kingdom | T: +44 (0)1423 817-733 | E: info@neoss.com