

M-VET THERAPY LASER

The Best Just Got Better

The M-VET Therapy Laser is a result of the scientific and technological know-how of over thirty years of research and development in a new generation of therapeutic devices entirely dedicated to animal health.

The benefits of the MLS® synchronized pulse pattern are now enhanced through an optical fibre technology, which allows conveying higher intensity energy packets inside the tissues, in total safety.

M-VET Treatment

- Exerts a significant inhibitory action on the production of pro-inflammatory molecules
- Produces short and long-term analgesic effects, allowing a significant and progressive function recovery
- Improves the quality of the newly formed tissue and accelerates healing times

Features

- Peak Power 270W
- Flexible and resistant fiber optic
- 10.1" high resolution LCD touchscreen display
- Pre-set protocols for pets, horses, reptiles, birds and exotic animals
- On-screen treatment guides
- Program sections for laser acupuncture and intraoral pathologies
- Customizable treatment parameters
- Variable treatment times
- Automatic calculation of the emitted energy
- Patient treatment log of up to 500 profiles with ability to export via USB port

Specifications

LASER TYPE	Class IV
WAVELENGTH	808nm and 905nm
M-VET DIMENSIONS	14 x 12 x 6 in, 9 lbs
TROLLEY DIMENSIONS	20 x 20 x 37 in, 42 lbs
POWER SUPPLY	100-240 V AC 50/60 Hz 51-68 VA



4 Interchangeable Lenses


The handpiece connection system ensures automatic adaptation of parameters.




DERMATOLOGY
Stimulation of tissue repair processes.



ORAL PATHOLOGIES
Analgesic, anti-inflammatory action and stimulation of tissue repair.



MUSCULOSKELETAL SYSTEM
Effective anti-inflammatory action through the reduction of inflammatory mediators.



LASERPUNCTURE
Stimulation of acupuncture points to rebalance the flow

Better Medicine. Better Business.

800.889.4184 | www.celasers.com

