

Ref: SB/4592


07 October 2020




Dear customer

Change Notification – TC78 / Change to 100ml Infusion bags

Following the change notification sent out in July 2020 with reference to printed 100ml infusion bags being temporarily supplied, the supply issue has been resolved and as such we will now be reverting back to unprinted 100ml infusion bags. Our supplier has made us aware that due to a new moulding tool there is a change in the position of the Luer lock and spike ports. These two ports have swapped position.



Current



New

The bags supplied are made of the same material and are from the same manufacturer.

This is a permanent change with immediate effect, with all orders now being supplied with the new bags. The following products will be impacted by this change:

| Product Description | Part No/Code |
|---------------------|--------------------------|
| TSB 100ml/100ml | 4.0194 & 4.0303 / 218010 |
| TSB 50ml/100ml | 4.0195 & 4.0304 / 218030 |
| D/S TSB 50ml/100ml | 4.0214 & 4.0319 / 218020 |
| D/S TSB 100ml/100ml | 4.0215 & 4.0320 / 218040 |
| IRR TSB 100ml/100ml | 5.0295 / 145000 |

If you require any further information relating to this notification please do not hesitate to contact the quality team via email: quality@cherwell-labs.co.uk

Yours sincerely

A handwritten signature in blue ink, appearing to read 'S Brimble'.

Steven Brimble
Quality Manager