

COVID-19 TESTING

Healthcare Providers

COMPREHENSIVE TESTING SOLUTION

Biodesix, a renowned diagnostic company with expertise in lung disease has recently launched a comprehensive, COVID-19 testing program for healthcare providers. We offer an end-to-end customized testing solution which includes specimen collection kits, coordinated logistics, billing support, and reliable EUA-FDA authorized tests in partnership with Bio-Rad Laboratories. Biodesix operates two high-complexity CAP/COLA accredited, CLIA and NYS CLEP certified laboratories with high capacity for COVID-19 testing.

WE ARE A TRUSTED PARTNER YOUR TEAM CAN RELY ON FOR COVID-19 TESTING

2 TESTS



FOR SARS-COV-2 DETECTION

24 HOURS

FROM SAMPLE RECEIPT TO RESULTS DELIVERY

22,000+

SPECIMENS TESTED IN FIRST
3 MONTHS POST-LAUNCH

COVID-19 TESTS PERFORMED BY BIODESIX	BIO-RAD SARS-COV-2 DDPCR TEST	PLATELIA SARS-COV-2 TOTAL AB TEST
What type of test is this?	MOLECULAR/DIAGNOSTIC	SEROLOGY/ANTIBODY
What does the test detect?	ACTIVE SARS-CoV-2 INFECTION	PREVIOUS SARS-CoV-2 INFECTION
What does the test measure?	SARS-COV-2 N GENE (N1 & N2) AND THE HUMAN RPP30 GENE	ANTIBODIES TO SARS-COV-2 (IGM, IGG & IGA)
What is the clinical use for the test?	IDENTIFY WHO SHOULD QUARANTINE	IDENTIFY WHO MAY HAVE PREVIOUSLY BEEN EXPOSED OR INFECTED
How reliable is the test data?	100% CONCORDANCE WITH CDC ASSAY symptomatic individuals	97.50% SENSITIVITY 99.56% SPECIFICITY >8 days post onset of symptoms
How are samples collected?	NASAL SWABS 	BLOOD DRAW 
What is the list price?	Medicare has set the price of \$100. The test list price is \$100.	Medicare has set the price of \$100. The test list price is \$100.

COVID-19 TESTING

End-to-End Customized Solution

We recognize everyone's COVID-19 situation is unique, which is why we wanted to provide you with a list of customizable steps in our COVID-19 testing program to fit your needs.



TEST ORDERING

- COVID-19 Tests** – option to order both the molecular and antibody test together or individually depending on your testing needs.



SPECIMEN COLLECTION



- Specimen Collection Kits** – we provide nasal swabs and blood collection kits for molecular and serology testing, respectively.
- Mobile Specimen Collection** – we can facilitate scheduling a healthcare provider to perform centralized specimen collection as needed.



SHIPPING



- Scheduled Courier Pick-Up (Colorado & Kansas)** – we can ensure same-day specimen pick-up with scheduled local couriers.
- FedEx Overnight Shipping (US)** – we use FedEx for overnight shipment of specimens to the Biodesix laboratory.
- Open All Week** – we accept specimens for testing 7-days per week to ensure we deliver results as soon as possible.



RESULT REPORTING

- Individual & Batch Test Reporting** – we offer two options for test result reporting that are HIPAA compliant.
- COVID-19 Demographic Report** – we can provide a demographic coversheet for a snapshot of the impact to your organization.



BILLING

- Biodesix is able to bill you, except as required by State or Federal regulations, or the patient's insurance for COVID-19 testing.