



Clinic RAMQ configuration form

Here are the information needed to set up a FMG clinic in the Omnimed EMR. **This information must be obtained from the RAMQ.** This interface will allow the client to receive a daily list of patients who are registered with the RAMQ and to consult the FMG or off FMG registration of a patient in the EMR as well as his vulnerability codes up to a maximum of three codes per patient.

| Full name of the site (FMG and off FMG) | FMG reference number (1XXX) |
|---|-----------------------------|
| | |

Name of your FMG site or off FMG and the access code AIP (ID MACHINE) of your site (AIPXXXXxx) as well as temporary password provided by the RAMQ. (e.g. nature2c)

| Name of the site | Access code (AIP) | Temporary password |
|------------------|-------------------|--------------------|
| | | |

FMG ONLY - Name of other FMG sites and respective AIP access codes (AIPXXXXxx) as well as temporary password provided by the RAMQ. (e.g. nature2c)

| Name of the site | Access code (AIP) | Temporary password |
|------------------|-------------------|--------------------|
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |

| FMG institution number (5XXXX) |
|--------------------------------|
| |

First and last name of a user with a valid CE code* and the password associated with the CE code to access the RAMQ. This password can be reset by the user on the RAMQ website as soon as the configuration is completed.

| Last and first name | CE Code | RAMQ Password |
|---------------------|---------|---------------|
| | | |

The CE code used is valid for the following doctors.

| Last and first name | CE Code |
|---------------------|---------|
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |

*If you need to apply for an AIP code from the RAMQ, please note that you only need one identifier.

The identified CE code must have access and permissions for the entire group of physicians. If the CE code doesn't have access to certain physicians, the synchronisation won't be authorized and won't work for the unauthorized physicians. Only one CE code is used for this synchronization step.